

IsoSep

product number

60/84

“T-antigen“-APE-HSA

Gal β 1-3GalNAc α -O-APE-HSA

The oligosaccharide is linked to HSA lysine residues via a spacer using the isothiocyanate coupling method. The spacer used is *p*-aminophenylethyl (APE). The spacer is glycosidically linked to the reducing monosaccharide unit of the oligosaccharide.

The conjugate normally contain 10-20 moles oligosaccharide per mole HSA.

The conjugate is supplied as freeze dried powder. Storage: -18° C.